

Exinpure DNA Spin Columns with Collection Tube

CAT # Exin-5005

Product Details:

Spin Column and collection tube (blue tube, clear ring, clear collection), Spin Columns are packed with silica membranes to rapidly purify nucleic acid using a silica-based extraction method. The columns can be used to purify single- or double-stranded DNA or RNA from contaminants such as salts, solvents, enzymes or proteins. Nucleic acids are bind to the silica resins under high chaotropic salt concentrations, salts and impurities are removed during the wash step. Pure DNA/RNA is eluted from the membrane using water or TE buffer. This isolation method of DNA/RNA is faster and easier to perform than other organic-based extraction methods. The EZ-10 Spin Column complete with attached lid and with 2-ml collection tubes, it is really convenient to use.

Product details at-a-glance / Technical Specifications

Total Product Size	100UNIT
Individual Container	Size:100UNIT
Shipping Conditions	RT

ORDERING INFORAMTION

CAT NO	PK / SIZE
Exin-5005	100 /pkt

